

We offer various packaging (protein concentration, activity, etc.) if necessary.

Data sheet

Enzyme	: Enolase YK1
--------	----------------------

Code ; ENL-97-01

Lot No. ;

Protein conc. ; mg/ml

Volume ; ml

Form ; 20 mM Tris-HCl (pH 7.5), 20 mM KCl


Storage ; -20 °C

Activity ; U/ml

Notes ; For research use only

Activity measurement

Reaction mix (20 mM Potassium phosphate buffer (pH 7.0), 1 mM 2-Phosphoglycerate, 0.3 mM NADH, 5 mM MgCl₂, 5 mM ADP, 10 U/ml Pyruvate kinase, 10 U/ml Lactate dehydrogenase and appropriate amount of the enzyme) was incubated at 50 °C and A₃₄₀ was monitored. One unit is defined as the amount of enzyme oxidizing 1 μmol of NADH ($\epsilon_{340}=6.22 \text{ mM}^{-1} \text{ cm}^{-1}$) per 1 minute using 2-Phosphoglycerate as a substrate.

pH profileThermostability